



Customized solutions for cell cultures.

Protocol

Production of **TCX6D Medium** from powder

Production of liquid medium with the TCX6D Powder Kit

Material:

- We recommend preparing the whole Powder Kit in a single batch! For that, please adjust the amounts/volumes per L given in this protocol according to your container/batch size!
- TCX6D Powder (27.01 g/L; Cat.No. 1070-XXXXDPM)
- ESK-Supplement (2.00 mL/L)
- approx. 1 L H₂O per L medium (WFI or equivalent quality)
- 7.50 mL/L 1 M NaOH Ph. Eur.
- 2.10 g/L NaHCO₃ Ph. Eur.
- 0.10 mL/L LONG® R3 IGF-I (Cat. No. 1006-XXXX);
Alternatively, 2 mL/L Growth hormone supplement (Cat. No. 1005-XXXX)
The values are corresponding to the use of 0.10 mg/L LONG® R3 IGF-I; alternatively: 10 mg/L rInsulin.
- We recommend wearing a dust mask during preparation!



Visual control:

- A. Container **Sealed and without any damage.**
- B. Appearance **Free flowing powder** (record color).




Check:















Color: _____

Procedure:

Check:

1.		Fill 0.8 L per 1 L (80% v/v) final medium 15-35°C water (WFI or equivalent quality) into the stirred tank/blending vessel.	<input type="checkbox"/>
2.		Start the stirrer of the system. Due to foam formation during medium production, the vortex should not reach the stirrer.	<input type="checkbox"/>
3.		Add 27.01 g/L of TCX6D powder slowly to the stirred water. Avoid clumping. Note: <i>We recommend preparing the whole powder kit at once.</i>	<input type="checkbox"/>

4.		Rinse the empty medium container with a suitable amount of water (WFI or equivalent quality) and pour liquid into the stirred tank.	<input type="radio"/>
5.		Stir for 30 minutes . Note: <i>The powder will not be completely dissolved at this stage.</i>	<input type="radio"/>
6.		Add 7.50 mL/L of 1 M NaOH to the stirred tank. <i>Adjust volume according to batch size.</i>	<input type="radio"/>
7.		Stir for 30 minutes . Note: <i>The powder should be completely dissolved and the solution should be clear.</i>	<input type="radio"/>
8.		Add 2.10 g/L NaHCO₃ to the stirred tank. <i>Adjust volume according to batch size.</i>	<input type="radio"/>
9.		Add an appropriate volume of water (WFI or equivalent quality) to reach the final volume. <i>Final volume depends on batch/container size!</i>	<input type="radio"/>
10.		Stir for ~ 10 minutes with lid closed. Note: <i>The solution should be clear, without precipitates. If not, stepwise increase mixing time by up to further 10 min.</i>	<input type="radio"/>
11.		Add 2.00 mL/L ESK Supplement to the stirred tank. <i>Adjust volume according to batch size.</i>	<input type="radio"/>

12.		<p>Add 0.10 mL/L LONG® R3 IGF-I (alternatively: 2.00 mL/L Growth hormone supplement)</p> <p>Note: <i>Adjustment of growth hormone concentration for optimization is possible, but depends on used cell line and application.</i></p>	<input type="radio"/>
13.		<p>Stir for ~ 10 minutes with lid closed.</p>	<input type="radio"/>
14.		<p>Check pH (pH 7.1 - pH 7.6) and osmolality (283 - 343 mOsmol/kg).</p>	<input type="radio"/> <input type="radio"/>
15.		<p>The Medium can now be sterile filtered (0.45 µm + 0.1µm) and bottled.</p>	<input type="radio"/>

Change History:

Revision	Date	Author	Comment/Description
01	31.03.2022	AWU	Initial version

www.xell.de
info@xell.de

Xell AG
Alte Verler Str. 1
33689 Bielefeld
Germany

For further information or assistance
contact us.

Fon: +49 (0)521 96989-200
Fax: +49 (0)521 96989-201

